MICROBIOLOGIA FUNCIONAL (2019-2020)

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
(<http://www.ruf.rice.edu/~bioslabs/methods/microscopy/cellcounting.html>)

To prepare the counting chamber the mirror-like polished surface is carefully cleaned with lens paper. The coverslip is also cleaned. Coverslips for counting chambers are specially made and are thicker than those for conventional microscopy, since they must be heavy enough to overcome the surface tension of a drop of liquid. The coverslip is placed over the counting surface prior to putting on the cell suspension. The suspension is introduced into one of the V-shaped wells with a pasteur or other type of pipet. The area under the coverslip fills by capillary action. Enough liquid should be introduced so that the mirrored surface is just covered. The charged counting chamber is then placed on the microscope stage and the counting grid is brought into focus at low power.

It is essential to be extremely careful with higher power objectives, since the counting chamber is much thicker than a conventional slide. The chamber or an objective lens may be damaged if the user is not careful. One entire grid on standard hemacytometers with Neubauer rulings can be seen at 40x (4x objective). The main divisions separate the grid into 9 large squares (like a tic-tac-toe grid). Each square has a surface area of one square mm, and the depth of the chamber is 0.1 mm. Thus the entire counting grid lies under a volume of 0.9 mm-cubed.

Suspensions should be dilute enough so that the cells or other particles do not overlap each other on the grid, and should be uniformly distributed. To perform the count, determine the magnification needed to recognize the desired cell type. Now systematically count the cells in selected squares so that the total count is 100 cells or so (number of cells needed for a statistically significant count). For large cells this may mean counting the four large corner squares and the middle one. For a dense suspension of small cells you may wish to count the cells in the four 1/25 sq. mm corners plus the middle square in the central square. Always decide on a specific counting pattern to avoid bias. For cells that overlap a ruling, count a cell as "in" if it overlaps the top or right ruling, and "out" if it overlaps the bottom or left ruling.

Here is a way to determine a particle count using a Neubauer hemocytometer. Suppose that you conduct a count as described above, and count 187 particles in the five small squares described. Each square has an area of 1/25 mm-squared (that is, 0.04 mm-squared) and depth of 0.1 mm. The total volume in each square is $(0.04) \times (0.1) = 0.004$ mm-cubed. You have five squares with combined volume of $5 \times (0.004) = 0.02$ mm-cubed. Thus you counted 187 particles in a volume of 0.02 mm-cubed, giving you $187 / (0.02) = 9350$ particles per mm-cubed. There are 1000 cubic millimeters in one cubic centimeter (same as a milliliter), so your particle count is 9,350,000 per ml.

Cells are often large enough to require counting over a larger surface area. For example, you might count the total number of cells in the four large corner squares plus the middle combined. Each square has surface area of 1 mm-squared and a depth of 0.1 mm, giving it a volume of 0.1 mm-cubed. Suppose that you counted 125 cells (total) in the five squares. You then have 125 cells per 0.5 mm-cubed, which is 250 cells/mm-cubed. Again, multiply by 1000 to determine cell count per ml (250,000).

Sometimes you will need to dilute a cell suspension to get the cell density low enough for counting. In that case you will need to multiply your final count by the dilution factor. For example, suppose that for counting you had to dilute a suspension of *Chlamydomonas* 10 fold. Suppose you obtained a final count of 250,000 cells/ml as described above. Then the count in the original (undiluted) suspension is $10 \times 250,000$ which is 2,500,000 cells/ml.

	Objetivos
--	------------------

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Use of Extracellular Acidification for the Rapid Testing of Ethanol Tolerance in Yeasts

J. Jiménez

Department of Genetics, Faculty of Biology, University of Sevilla,
Apto 1095, Sevilla, Spain

N. van Uden*

Laboratory of Microbiology, Gulbenkian Institute of Science, 2781 Oeiras
Codex, Portugal

Accepted for publication February 12, 1985

The overall kinetics of the inhibition by ethanol of growth and fermentation in *Saccharomyces cerevisiae* and other yeasts is composite and depends on the contribution of a number of underlying inhibitory mechanisms (for a review, see van Uden¹). Among these, the interference of ethanol with membrane-bound processes plays a dominant role, including the exponential enhancement of thermal death² and of passive proton influx,³ and the exponential noncompetitive inhibition of the transport of glucose,⁴ maltose,⁵ ammonium,⁶ and amino acids.⁷ The exponential constant of ethanol-enhanced passive proton diffusion across the plasma membrane of yeast can be estimated rapidly and easily. It occurred to us that this might constitute a rapid method for the preliminary evaluation of ethanol tolerance in yeasts, such as would be convenient to use when large numbers of mutant and other strains have to be screened. To test its feasibility we verified, using a number of strains, whether a correlation existed between the values of the exponential constant (k) of ethanol-enhanced proton influx and the ethanol concentration (K_i) that reduced the specific growth rate by 50%.

MATERIAL AND METHODS

Information on the 11 yeast strains used in the study is listed in Table I. The strains were grown with magnetic stirring at 25°C in thermostated Erlenmeyer flasks containing a liquid medium of the following composition: glucose, 2%; yeast extract, 0.5%; peptone, 1%; with or without added ethanol.

The specific growth rates were calculated in the usual way from the linear parts of semilog plots of the optical densities of samples against the sampling times. To estimate K_i , the concentration of ethanol

that reduced the specific growth rate by 50%, growth experiments were performed using batteries of flasks containing the above medium to which ethanol had been added to final concentrations of 0%, 1%, 2%, 3%, 4%, 5%, 6%, 7%, 8%, 9%, and 10% (v/v). Estimates of K_i were obtained by visual inspection of plots of the specific growth rates against ethanol concentrations. The K_i values finally retained were averages of from three to five sets of experiments done with each strain (Table I).

For estimating k , the exponential constant of ethanol-enhanced proton diffusion, the yeasts were grown in 500 ml medium into early stationary phase, centrifuged, washed twice with distilled water, and resuspended in 10 ml distilled water. Of this suspension, 0.5 ml was placed in a 20-ml glass vial together with 0.5 ml of a 20% glucose solution and distilled water with or without ethanol at the desired concentration to a total volume of 4.5 ml. The ethanol concentrations used were 0%, 2%, 4%, 6%, 8%, and 10% (v/v). The vials were screw-capped and incubated at 25°C with mechanical shaking during 75 min, after which the pH was measured with a standard pH meter PHM62 (Radiometer, Copenhagen).

In preliminary experiments done with all strains using the pH meter connected to a flat-bed recorder LKB 6500, it was established that after 75 min the final pH values were stable while glucose was still present ensuring active proton efflux (without which the pH would start to rise again).

RESULTS AND DISCUSSION

For determining the exponential enhancement constants of proton influx, the following relations were used.³ When yeast cells are suspended in water containing glucose, protons are extruded and the extracellular pH drops to a final value following a char-

* To whom all correspondence should be addressed.

Table I. Yeast strains: identity, origin, and parameters of ethanol toxicity.

Serial number	Strain	Species	Origin	K_i^* (M)	k^\dagger (M ⁻¹)
1	IFI 256	<i>Saccharomyces cerevisiae</i>	Wine from Galicia, Spain	1.24 (0.08)	0.55 (0.998)
2	IFI 82	<i>Saccharomyces cerevisiae</i>	Wine from Andalusia, Spain	1.13 (0.10)	0.80 (0.951)
3	ACA 4	<i>Torulaspota delbrueckii</i>	Grape must from Extremadura, Spain	1.22 (0.07)	0.83 (0.999)
4	ACA 174	<i>Saccharomyces cerevisiae</i>	Grape must from Andalusia, Spain	1.37 (0.09)	0.26 (0.967)
5	RBV-DSA	<i>Saccharomyces cerevisiae</i>	Meiotic product of a distillery strain	1.26 (0.11)	0.59 (0.995)
6	K5-SA	<i>Saccharomyces cerevisiae</i>	Genetic strain	0.96 (0.10)	1.01 (0.990)
7	FIF 414	<i>Torulaspota pretoriensis</i>	Flor of Sherry wine from Andalusia, Spain	1.50 (0.08)	0.15 (0.862)
8	X30-36	<i>Saccharomyces cerevisiae</i>	Genetic strain	0.96 (0.16)	0.70 (0.963)
9	FSP 6	<i>Torulaspota pretoriensis</i>	Meiotic product of strain 7	1.37 (0.06)	0.45 (0.976)
10	FXH 1	hybrid	8 × 9	1.16 (0.15)	0.51 (0.921)
11	IGC 3507	<i>Saccharomyces cerevisiae</i>	Respiration-deficient mutant	1.26	0.55 (0.991)

* This is the concentration of added ethanol that reduced the specific growth rate by 50% (standard deviation in parentheses).

† This is the exponential constant of ethanol-enhanced proton diffusion across plasma membrane (correlation coefficient from linear regression).

acteristic acidification curve. At the final pH value the active extrusion of protons is balanced by passive proton influx. Ethanol exponentially enhances the diffusion constant of the proton influx:

$$C_X = C_0 e^{kX} \quad (1)$$

where C_X and C_0 are the values of the diffusion constants at concentration X and 0 of ethanol, respectively, and k is the exponential enhancement constant. Consequently, ethanol increases the final pH value according to

$$\text{pH}_f = \log_{10} C_0 - \log_{10} k_p + kX \log_{10} e \quad (2)$$

where k_p is the rate of active proton extrusion. Under the condition that $(\log_{10} C_0 - \log_{10} k_p)$ is a constant or nearly so over the range of experimental pH_f values, eq. (2), if valid, should be linear and estimates of k may be obtained from the slope (for a more detailed theoretical treatment see Leão and van Uden³).

In 10 of the 11 strains the relations between final pH and ethanol concentrations were indeed linear with correlation coefficients ranging from 0.951 to 0.999, while in strain 7 the linear fit was poorer (correlation coefficient 0.862). Using eq. (2), estimates for k were calculated from the slopes of the regression lines (Table I).

As is shown in Fig. 1, the relationship in the 11 strains between K_i (concentration of added ethanol that reduces the specific growth rate by 50%) and k (exponential constant of ethanol-enhanced proton dif-

fusion) was roughly linear (correlation coefficient of 0.848).

We conclude that it is worthwhile to determine k with the rapid method described above when large numbers of strains have to be screened and that the values obtained provide preliminary information on the probable effects of ethanol on growth. Whether a similar correlation applies to other toxic effects of ethanol on yeasts remains to be seen.


Figure 1. Relations between K_i (concentration of added ethanol that reduced the specific growth rate of yeast growth by 50%) and k (the exponential constant of ethanol-enhanced proton diffusion across the yeast plasma membrane). Numbers refer to yeast strains as listed in Table I.

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